










Prepare patient for collection

<p>1</p> <p>Confirm that the patient has not opened the Quantisal™ packaging prior to initiation of preparation for collection.</p> 	<p>2</p> <p>To ensure that the patient's mouth is clear, instruct him/her to swish and spit with water at least 10 minutes prior to collection. The patient is to refrain from consumption of food (including gum, candy, etc.) or beverage (including water, coffee, soda, etc.) during the 10 minutes prior to collection.</p> 	<p>3</p> <p>Ask the patient to open the Quantisal™ package, and remove the collection device and transport tube from the package.</p> <p>Confirm the expiration date has not passed on the Quantisal™ packaging</p> 
<p>4</p> <p>Instruct the patient to wash hands thoroughly with soap and water. Request that they wet their hands first, apply soap, and rub hands vigorously for at least 15 seconds. Then, request that the patient rinse hands with water and use disposable towel to dry. Towel should be used to turn off faucet.</p> 	<p>5</p> <p>Have the patient open the pouch. Request that the patient provide visual evidence that Volume Adequacy Indicator IS NOT blue after collection device is removed from pouch. Patient should not set the collection device down once it has been removed from packaging.</p> 	<p>6</p> <p>Request that the patient accumulate saliva (refrain from swallowing) in his/her mouth before starting collection.</p> 

Sample Collection


<p>1</p> <p>Instruct the patient to position collection device under the tongue and close their mouth as they would with a thermometer.</p> 	<p>2 Ensure patient remains within line of sight during entire collection period.</p> <p>IMPORTANT: Advise the patient to avoid chewing or sucking the pad. The collection device should remain under the tongue until the Volume Adequacy Indicator turns completely blue.</p> <p>Please instruct the patient to not remove the collection device at the first indication of blue in the indicator window (see below). The device can be removed as soon as the patient provides visual evidence that the indicator has turned completely blue.</p> <p>If the indicator has not turned blue in 10 minutes, ask the patient to remove the device and discard. If available, the patient may attempt to collect a sample using a new device after drinking water and waiting 10 minutes prior to recollection.</p> <div style="display: flex; align-items: center; justify-content: center;"> <div style="text-align: center;"> <p>FAILED ATTEMPT</p>  </div> <div style="margin: 0 20px;">→</div> <div style="text-align: center;"> <p>SUCCESSFUL ATTEMPT</p>  </div> </div>
--	--

Preparation of sample prior to shipment

1


While still holding the collection device in one hand, instruct patient to pick up the transport tube in opposite hand and uncap by pushing up.

Do not spill or empty blue buffer liquid from the tube. Visually inspect transport tube to ensure that the blue buffer liquid is present.



2


Once transport tube is uncapped, instruct patient to insert the collection device into the uncapped transport tube and carefully replace the cap. The pad should be placed in the blue buffer liquid. Patient should not ingest blue buffer solution or place device in mouth after it has made contact with the buffer solution.



3

Instruct the patient to carefully place cap over the top of the collection device in the transport tube.

A “snap” should be heard once the cap is successfully placed on the transport tube. Once the cap has been placed on the transport tube, request that the patient mix gently by shaking the tube.




4

Request that the patient retrieve the seal containing a barcode from the package and ask that they write their name and date of birth on the seal.



5

Ask the patient to provide you with the barcode number (i.e. 81410349 as shown in the example below).



Request that the patient place the seal with name and date of birth included on the side of the transport tube.

6

Instruct patient to place the capped transport tube in the clear Aegis specimen bag and seal the bag closed.









7

Ask the patient to place the specimen bag in the US Postal Service (USPS) shipping package, seal the package closed, and place the USPS label on the outside of the package.



Instruct the patient to place the sealed USPS package in the mailbox and lift the red flag on the mailbox if appropriate.

IMPORTANT - BEFORE PROCEEDING - CHECK THE PATIENT'S IDENTIFICATION.

Prepare patient for collection

<p>1</p> <p>Confirm that the patient has not opened the Quantisal™ packaging prior to initiation of preparation for collection.</p> 	<p>2</p> <p>Verify that the patient's mouth is clear and that they have refrained from consumption of food (including gum, candy, etc.) or beverage (including water, coffee, soda, etc.) 10 minutes prior to collection.</p> 	<p>3</p> <p>Ask the patient to open the Quantisal™ package, and remove the collection device and transport tube from the package.</p> <p>Confirm the expiration date has not passed on the Quantisal™ packaging</p> 
<p>4</p> <p>Instruct the patient to wash hands thoroughly with soap and water. Request that they wet their hands first, apply soap, and rub hands vigorously for at least 15 seconds. Then, request that the patient rinse hands with water and use disposable towel to dry. Towel should be used to turn off faucet.</p> 	<p>5</p> <p>Have the patient open the pouch. Request that the patient provide visual evidence that Volume Adequacy Indicator IS NOT blue after collection device is removed from pouch. Patient should not set the collection device down once it has been removed from packaging.</p> 	<p>6</p> <p>Request that the patient accumulate saliva (refrain from swallowing) in his/her mouth before starting collection.</p> 

Sample Collection

<p>1</p> <p>Instruct the patient to position collection device under the tongue and close their mouth as they would with a thermometer.</p> 	<p>2 Ensure patient remains within line of sight during entire collection period.</p> <p>IMPORTANT: Advise the patient to avoid chewing or sucking the pad. The collection device should remain under the tongue until the Volume Adequacy Indicator turns completely blue.</p> <p>Please instruct the patient to not remove the collection device at the first indication of blue in the indicator window (see below). The device can be removed as soon as the patient provides visual evidence that the indicator has turned completely blue.</p> <p>If the indicator has not turned blue in 10 minutes, ask the patient to remove the device and discard. If available, the patient may attempt to collect a sample using a new device after drinking water and waiting 10 minutes prior to recollection.</p> 
--	--


Collector Instructions

Preparation of sample prior to shipment

1


While still holding the collection device in one hand, instruct patient to hold the transport tube and uncap by pushing up.

Do not spill or empty blue buffer liquid from the tube. Visually inspect transport tube to ensure that the blue buffer liquid is present.



2


Once transport tube is uncapped, instruct patient to insert the collection device into the uncapped transport tube and carefully replace the cap. The pad should be placed in the blue buffer liquid. Patient should not ingest blue buffer solution or place device in mouth after it has made contact with the buffer solution.



3


Instruct the patient to carefully place cap over the top of the collection device in the transport tube.

A “snap” should be heard once the cap is successfully placed on the transport tube. Once the cap has been placed on the transport tube, request that the patient mix gently by shaking the tube.




4

Request that the patient retrieve the seal containing a barcode from the package and ask that they write their name and date of birth on the seal.



5

Request that the patient place the seal with name and date of birth included on the side of the transport tube.



6

Instruct patient to place the capped transport tube in the clear Aegis specimen bag and seal the bag closed.



7

Ask the patient to place the specimen bag in the US Postal Service (USPS) shipping package, seal the package closed, and place the USPS label on the outside of the package.

Instruct the patient to place the sealed USPS package in the mailbox and lift the red flag on the mailbox if appropriate.

Important Instructions

PLEASE READ THIS FIRST

This package contains everything needed for an oral fluid collection: collection device, barcode label, and USPS mailer.

This test has been ordered by your treatment provider and is part of your treatment program.

It is very important that you read and follow the steps outlined below.

Prior to Your Test

Please do not open the oral swab package until your appointment time and are instructed to do so as opening the swabs could invalidate your test.

Once you are notified of your telemedicine appointment time, please follow these steps to ensure a smooth collection:

- Drink some water 60 minutes prior to your appointment time.
- However, please do not eat any food (including gum, candy, etc.) or drink any beverages (including water, coffee, soda, etc.) 10 minutes prior to your appointment time.
- Please minimize distractions during your appointment time as your full attention will be appreciated.
- Please make sure you have your identification and the contents of the testing packet in front of you for your appointment. Please do not open the oral swab package until asked to do so by your provider or person witnessing the collection.
- Your provider or person witnessing the collection will walk you through the step by step instructions to complete the oral fluid collection.
- Once the collection is complete and the package is ready to be shipped, please place in your mailbox and lift the red flag to signal to the mail courier a package is ready for delivery.
- If you have any questions and wish to speak with an Aegis representative, please call 1-800-533-7052.